
thermo scientific

EliA™

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
- "Screening for CTD-associated antibodies by automated immunoassay", Willems et al. 2018⁹
- "The association of solid-phase assays to immunofluorescence increases the diagnostic accuracy for ANA screening in patients with autoimmune rheumatic diseases", Bizzaro et al. 2018⁴

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- “A comparison of a fluorescence enzyme immunoassay versus IIF for initial screening of CTD”, Orme et al. 2018¹⁰
- “Measurement of antinuclear antibodies and their fine specificities: time for a change in strategy?”, Otten et al. 2017⁵
- “Comparison of the clinical utility of the EIA CTD Screen to IIF on Hep-2 cells”, Robier et al. 2016⁷

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- "2019 EULAR/ACR classification criteria for systemic lupus erythematosus", Aringer et al. 2019³
- "EULAR/ACR SLE classification criteria item performance", Aringer et al. 2021¹⁴

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
- "Sydney Criteria" as part of the "11th International Congress on antiphospholipid antibodies", Miyakis et al. 2006 ²¹
- "International Consensus Guidelines on Anticardiolipin and Anti-β2-Glycoprotein I Testing", 2012 ²²

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
- "2010 ACR/EULAR classification criteria for rheumatoid arthritis", Aletaha et al. 2010²⁶
- "Relationship between rheumatoid factor isotypes and IgG anti-cyclic citrullinated peptide antibodies", Jaskowski et al. 2010²⁷
- "Determination of Autoantibody Isotypes Increases the Sensitivity of Serodiagnostics in Rheumatoid Arthritis", Sieghart et al. 2018²⁸

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- "Revised 2017 international consensus on testing of ANCAs in GPA and MPA", 2017⁴³

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
- "ECCO-ESGAR Guideline for Diagnostic Assessment in IBD Part 1: Initial diagnosis, monitoring of known IBD, detection of complications", Maaser et al. 2019⁵²
- "World Gastroenterology Organisation Global Guidelines Inflammatory Bowel Disease: Update August 2015", Bernstein et al. 2016⁵³

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- “Guidelines for Diagnosing Coeliac Disease” from the European Society Paediatric Gastroenterology, Hepatology and Nutrition (ESPGHAN), Husby et al. 2020⁶³
- “European Society for the Study of Coeliac Disease (ESsCD) guideline for coeliac disease and other gluten-related disorders”, Al-Toma et al. 2019⁶⁴
- “Accuracy of a no-biopsy approach for the diagnosis of coeliac disease across different adult cohorts”, Penny et al. 2020⁶⁵

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- "EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis", European Association for the Study of the Liver, 2017⁶⁹
- "EASL Clinical Practice Guidelines: Autoimmune hepatitis", European Association for the Study of the Liver, 2015⁷⁰
- "Diagnosis and Management of Autoimmune Hepatitis: Current Status and Future Directions", Czaja, 2016⁷¹

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- "European Thyroid Association Guideline for the Management of Graves' Hyperthyroidism", Kahaly et al. 2018⁷⁴
- "European Thyroid Association Guidelines on the Diagnosis and Management of Central Hypothyroidism", Persani et al. 2018⁷⁵

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1. Gaubitz M. Epidemiology of connective tissue disorders. *Rheumatology (Oxford)*. 2006;45 Suppl 3(3):iii3-4.
2. Damoiseaux J, Olschowka N, Shoenfeld Y. EASI - European Autoimmunity Standardisation Initiative: facing the challenges of diagnostics in autoimmunity. *Clin Chem Lab Med*. 2018;56(10):1620-3.
3. Aringer M, Costenbader K, Daikh D, Brinks R, Mosca M, Ramsey-Goldman R, et al. 2019 European League Against Rheumatism/American College of Rheumatology classification criteria for systemic lupus erythematosus. *Ann Rheum Dis*. 2019;78(9):1151-9.
4. Bizzaro N, Brusca I, Previtali G, Alessio MG, Daves M, Platzgummer S, et al. The association of solid-phase assays to immunofluorescence increases the diagnostic accuracy for ANA screening in patients with autoimmune rheumatic diseases. *Autoimmun Rev*. 2018;17(6):541-7.
5. Otten HG, Brummelhuis WJ, Fritsch-Stork R, Leavis HL, Wisse BW, van Laar JM, et al. Measurement of antinuclear antibodies and their fine specificities: time for a change in strategy? *Clinical and Experimental Rheumatology*. 2017;35(3):462-70.
6. Robier C, Amouzadeh-Ghadikolai O. The clinical significance of borderline results of the Elia CTD Screen assay. *Clin Chem Lab Med*. 2018;56(12):2088-92.
7. Robier C, Amouzadeh-Ghadikolai O, Stettin M, Reich G. Comparison of the clinical utility of the Elia CTD Screen to indirect immunofluorescence on Hep-2 cells. *Clin Chem Lab Med*. 2016;54(8):1365-70.
8. van der Pol P, Bakker-Jonges LE, Kuijpers J, Schreurs MWJ. Analytical and clinical comparison of two fully automated immunoassay systems for the detection of autoantibodies to extractable nuclear antigens. *Clin Chim Acta*. 2018;476:154-9.
9. Willems P, De Langhe E, Claessens J, Westhovens R, Van Hoeyveld E, Poesen K, et al. Screening for connective tissue disease-associated antibodies by automated immunoassay. *Clin Chem Lab Med*. 2018;56(6):909-18.
10. Orme ME, Andalucia C, Sjolander S, Bossuyt X. A comparison of a fluorescence enzyme immunoassay versus indirect immunofluorescence for initial screening of connective tissue diseases: Systematic literature review and meta-analysis of diagnostic test accuracy studies. *Best Pract Res Clin Rheumatol*. 2018;32(4):521-34.
11. Op De Beeck K, Vermeersch P, Verschueren P, Westhovens R, Marien G, Blockmans D, et al. Detection of antinuclear antibodies by indirect immunofluorescence and by solid phase assay. *Autoimmun Rev*. 2011;10(12):801-8.
12. Tan EM, Feltkamp TE, Smolen JS, Butcher B, Dawkins R, Fritzler MJ, et al. Range of antinuclear antibodies in "healthy" individuals. *Arthritis Rheum*. 1997;40(9):1601-11.
13. Fritzler MJ. The antinuclear antibody test: last or lasting gasp? *Arthritis Rheum*. 2011;63(1):19-22.
14. Aringer M, Brinks R, Dorner T, Daikh D, Mosca M, Ramsey-Goldman R, et al. European League Against Rheumatism (EULAR)/American College of Rheumatology (ACR) SLE classification criteria item performance. *Ann Rheum Dis*. 2021.
15. Choi MY, Clarke AE, St Pierre Y, Hanly JG, Urowitz MB, Romero-Diaz J, et al. Antinuclear Antibody-Negative Systemic Lupus Erythematosus in an International Inception Cohort. *Arthritis Care Res (Hoboken)*. 2019;71(7):893-902.
16. Bossuyt X, Luyckx A. Antibodies to extractable nuclear antigens in antinuclear antibody-negative samples. *Clin Chem*. 2005;51(12):2426-7.
17. Pourmand N, Blomberg S, Ronnblom L, Karlsson-Parra A, Pettersson I, Wahren-Herlenius M. Ro 52kD autoantibodies are detected in a subset of ANA-negative sera. *Scand J Rheumatol*. 2000;29(2):116-23.
18. Shiboski CH, Shiboski SC, Seror R, Criswell LA, Labetoulle M, Lietman TM, et al. 2016 American College of Rheumatology/European League Against Rheumatism Classification Criteria for Primary Sjogren's Syndrome: A Consensus and Data-Driven Methodology Involving Three International Patient Cohorts. *Arthritis Rheumatol*. 2017;69(1):35-45.
19. Chan EK, Damoiseaux J, Carballo OG, Conrad K, de Melo Cruvinel W, Francescantonio PL, et al. Report of the First International Consensus on Standardized Nomenclature of Antinuclear Antibody HEp-2 Cell Patterns 2014-2015. *Front Immunol*. 2015;6:412.
20. Conrad K, Schöbler W, Hiepe F, Fritzler MJ. Autoantibodies in Systemic Autoimmune Diseases – A Diagnostic Reference. 3rd ed: Pabst Science Publishers; 2015.
21. Miyakis S, Lockshin MD, Atsumi T, Branch DW, Brey RL, Cervera R, et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost*. 2006;4(2):295-306.

22. Lakos G, Favaloro EJ, Harris EN, Meroni PL, Tincani A, Wong RC, et al. International consensus guidelines on anticardiolipin and anti-beta2-glycoprotein I testing: report from the 13th International Congress on Antiphospholipid Antibodies. *Arthritis Rheum.* 2012;64(1):1-10.
23. Saigal R, Kansal A, Mittal M, Singh Y, Ram H. Antiphospholipid antibody syndrome. *J Assoc Physicians India.* 2010;58:176-84.
24. Sciascia S, Cosseddu D, Montaruli B, Kuzenko A, Bertero MT. Risk Scale for the diagnosis of antiphospholipid syndrome. *Ann Rheum Dis.* 2011;70(8):1517-8.
25. Pengo V, Ruffatti A, Del Ross T, Tonello M, Cuffaro S, Hoxha A, et al. Confirmation of initial antiphospholipid antibody positivity depends on the antiphospholipid antibody profile. *J Thromb Haemost.* 2013;11(8):1527-31.
26. Aletaha D, Neogi T, Silman AJ, Funovits J, Felson DT, Bingham CO, 3rd, et al. 2010 rheumatoid arthritis classification criteria: an American College of Rheumatology/European League Against Rheumatism collaborative initiative. *Ann Rheum Dis.* 2010;69(9):1580-8.
27. Jaskowski TD, Hill HR, Russo KL, Lakos G, Szekanecz Z, Teodorescu M. Relationship between rheumatoid factor isotypes and IgG anti-cyclic citrullinated peptide antibodies. *J Rheumatol.* 2010;37(8):1582-8.
28. Sieghart D, Platzer A, Studenic P, Alasti F, Grundhuber M, Swiniarski S, et al. Determination of Autoantibody Isotypes Increases the Sensitivity of Serodiagnostics in Rheumatoid Arthritis. *Front Immunol.* 2018;9:876.
29. Smolen JS, Landewe R, Bijlsma J, Burmester G, Chatzidionysiou K, Dougados M, et al. EULAR recommendations for the management of rheumatoid arthritis with synthetic and biological disease-modifying antirheumatic drugs: 2016 update. *Ann Rheum Dis.* 2017;76(6):960-77.
30. Suresh E. Diagnosis of early rheumatoid arthritis: what the non-specialist needs to know. *Journal of the Royal Society of Medicine.* 2004;97(9):421-4.
31. Littlejohn EA, Monrad SU. Early Diagnosis and Treatment of Rheumatoid Arthritis. *Prim Care.* 2018;45(2):237-55.
32. Bizzaro N. Antibodies to citrullinated peptides: a significant step forward in the early diagnosis of rheumatoid arthritis. *Clin Chem Lab Med.* 2007;45(2):150-7.
33. Berglin E, Padyukov L, Sundin U, Hallmans G, Stenlund H, Van Venrooij WJ, et al. A combination of autoantibodies to cyclic citrullinated peptide (CCP) and HLA-DRB1 locus antigens is strongly associated with future onset of rheumatoid arthritis. *Arthritis Res Ther.* 2004;6(4):R303-8.
34. Schellekens GA, Visser H, de Jong BA, van den Hoogen FH, Hazes JM, Breedveld FC, et al. The diagnostic properties of rheumatoid arthritis antibodies recognizing a cyclic citrullinated peptide. *Arthritis Rheum.* 2000;43(1):155-63.
35. de Brito Rocha S, Baldo DC, Andrade LEC. Clinical and pathophysiologic relevance of autoantibodies in rheumatoid arthritis. *Adv Rheumatol.* 2019;59(1):2.
36. Jonsson T, Thorsteinnsson J, Valdimarsson H. Elevation of only one rheumatoid factor isotype is not associated with increased prevalence of rheumatoid arthritis--a population based study. *Scand J Rheumatol.* 2000;29(3):190-1.
37. Bas S, Perneger TV, Kunzle E, Vischer TL. Comparative study of different enzyme immunoassays for measurement of IgM and IgA rheumatoid factors. *Ann Rheum Dis.* 2002;61(6):505-10.
38. Bobbio-Pallavicini F, Caporali R, Alpini C, Avalle S, Epis OM, Klersy C, et al. High IgA rheumatoid factor levels are associated with poor clinical response to tumour necrosis factor alpha inhibitors in rheumatoid arthritis. *Ann Rheum Dis.* 2007;66(3):302-7.
39. Shakiba Y, Koopah S, Jamshidi AR, Amirzargar AA, Massoud A, Kiani A, et al. Anti-cyclic citrullinated peptide antibody and rheumatoid factor isotypes in Iranian patients with rheumatoid arthritis: evaluation of clinical value and association with disease activity. *Iran J Allergy Asthma Immunol.* 2014;13(3):147-56.
40. de Angelis V, Meroni PL. Rheumatoid Factors. In: Shoenfeld Y, Gershwin ME, Meroni PL, editors. *Autoantibodies (Second Edition)*: Elsevier; 2007. p. 755-62.
41. Mikuls TR, O'Dell JR, Stoner JA, Parrish LA, Arend WP, Norris JM, et al. Association of rheumatoid arthritis treatment response and disease duration with declines in serum levels of IgM rheumatoid factor and anti-cyclic citrullinated peptide antibody. *Arthritis Rheum.* 2004;50(12):3776-82.
42. Song YW, Kang EH. Autoantibodies in rheumatoid arthritis: rheumatoid factors and anticitrullinated protein antibodies. *QJM.* 2010;103(3):139-46.

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43. Bossuyt X, Cohen Tervaert JW, Arimura Y, Blockmans D, Flores-Suarez LF, Guillevin L, et al. Position paper: Revised 2017 international consensus on testing of ANCAs in granulomatosis with polyangiitis and microscopic polyangiitis. *Nat Rev Rheumatol*. 2017;13(11):683-92.
44. Csernok E, Moosig F. Current and emerging techniques for ANCA detection in vasculitis. *Nat Rev Rheumatol*. 2014;10(8):494-501.
45. Savage J, Gillis D, Benson E, Davies D, Esnault V, Falk RJ, et al. International Consensus Statement on Testing and Reporting of Antineutrophil Cytoplasmic Antibodies (ANCA). *Am J Clin Pathol*. 1999;111(4):507-13.
46. Toh B-H. Autoimmune Gastritis and Pernicious Anemia. In: Rose NR, Mackay IR, editors. *The Autoimmune Diseases (Sixth Edition)*: Academic Press; 2019. p. 833-47.
47. Devalia V, Hamilton MS, Molloy AM, British Committee for Standards in H. Guidelines for the diagnosis and treatment of cobalamin and folate disorders. *Br J Haematol*. 2014;166(4):496-513.
48. Toh BH, Alderuccio F. Pernicious anaemia. *Autoimmunity*. 2004;37(4):357-61.
49. Andres E, Serraj K. Optimal management of pernicious anemia. *J Blood Med*. 2012;3:97-103.
50. Varbanova M, Frauenschlager K, Malfertheiner P. Chronic gastritis - an update. *Best Pract Res Clin Gastroenterol*. 2014;28(6):1031-42.
51. Rusak E, Chobot A, Krzywicka A, Wenzlau J. Anti-parietal cell antibodies - diagnostic significance. *Adv Med Sci*. 2016;61(2):175-9.
52. Maaser C, Sturm A, Vavricka SR, Kucharzik T, Fiorino G, Annese V, et al. ECCO-ESGAR Guideline for Diagnostic Assessment in IBD Part 1: Initial diagnosis, monitoring of known IBD, detection of complications. *J Crohns Colitis*. 2019;13(2):144-64.
53. Bernstein CN, Eliakim A, Fedail S, Fried M, Geary R, Goh KL, et al. World Gastroenterology Organisation Global Guidelines Inflammatory Bowel Disease: Update August 2015. *J Clin Gastroenterol*. 2016;50(10):803-18.
54. Molodecky NA, Soon IS, Rabi DM, Ghali WA, Ferris M, Chernoff G, et al. Increasing incidence and prevalence of the inflammatory bowel diseases with time, based on systematic review. *Gastroenterology*. 2012;142(1):46-54.
55. Harkness EF, Harrington V, Hinder S, O'Brien SJ, Thompson DG, Beech P, et al. GP perspectives of irritable bowel syndrome--an accepted illness, but management deviates from guidelines: a qualitative study. *BMC Fam Pract*. 2013;14:92.
56. Tibble JA, Sigthorsson G, Foster R, Forgacs I, Bjarnason I. Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease. *Gastroenterology*. 2002;123(2):450-60.
57. van Rhee PF, Van de Vijver E, Fidler V. Faecal calprotectin for screening of patients with suspected inflammatory bowel disease: diagnostic meta-analysis. *BMJ*. 2010;341:c3369.
58. Chew TS, Mansfield JC. Can faecal calprotectin predict relapse in inflammatory bowel disease: a mini review. *Frontline Gastroenterol*. 2018;9(1):23-8.
59. Costa F, Mumolo MG, Ceccarelli L, Bellini M, Romano MR, Sterpi C, et al. Calprotectin is a stronger predictive marker of relapse in ulcerative colitis than in Crohn's disease. *Gut*. 2005;54(3):364-8.
60. Roseth AG, Aadland E, Grzyb K. Normalization of faecal calprotectin: a predictor of mucosal healing in patients with inflammatory bowel disease. *Scand J Gastroenterol*. 2004;39(10):1017-20.
61. Schoepfer AM, Beglinger C, Straumann A, Trummel M, Renzulli P, Seibold F. Ulcerative colitis: correlation of the Rachmilewitz endoscopic activity index with fecal calprotectin, clinical activity, C-reactive protein, and blood leukocytes. *Inflamm Bowel Dis*. 2009;15(12):1851-8.
62. Theede K, Holck S, Ibsen P, Kallemose T, Nordgaard-Lassen I, Nielsen AM. Faecal Calprotectin Predicts Relapse and Histological Mucosal Healing in Ulcerative Colitis. *Inflamm Bowel Dis*. 2016;22(5):1042-8.
63. Husby S, Koletzko S, Korponay-Szabo I, Kurppa K, Mearin ML, Ribes-Koninckx C, et al. European Society Paediatric Gastroenterology, Hepatology and Nutrition Guidelines for Diagnosing Coeliac Disease 2020. *J Pediatr Gastroenterol Nutr*. 2020;70(1):141-56.
64. Al-Toma A, Volta U, Auricchio R, Castillejo G, Sanders DS, Cellier C, et al. European Society for the Study of Coeliac Disease (ESsCD) guideline for coeliac disease and other gluten-related disorders. *United European Gastroenterol J*. 2019;7(5):583-613.

65. Penny HA, Raju SA, Lau MS, Marks LJ, Baggus EM, Bai JC, et al. Accuracy of a no-biopsy approach for the diagnosis of coeliac disease across different adult cohorts. *Gut*. 2020.
66. Choung RS, Khaleghi S, Cartee AK, Marietta EV, Larson JJ, King KS, et al. Community-Based Study of Celiac Disease Autoimmunity Progression in Adults. *Gastroenterology*. 2020;158(1):151-9 e3.
67. Rubio-Tapia A, Hill ID, Kelly CP, Calderwood AH, Murray JA, American College of G. ACG clinical guidelines: diagnosis and management of celiac disease. *Am J Gastroenterol*. 2013;108(5):656-76; quiz 77.
68. Werkstetter KJ, Korponay-Szabo IR, Popp A, Villanacci V, Salemme M, Heilig G, et al. Accuracy in Diagnosis of Celiac Disease Without Biopsies in Clinical Practice. *Gastroenterology*. 2017;153(4):924-35.
69. EASL. EASL Clinical Practice Guidelines: The diagnosis and management of patients with primary biliary cholangitis. *J Hepatol*. 2017;67(1):145-72.
70. European Association for the Study of the L. EASL Clinical Practice Guidelines: Autoimmune hepatitis. *J Hepatol*. 2015;63(4):971-1004.
71. Czaja AJ. Diagnosis and Management of Autoimmune Hepatitis: Current Status and Future Directions. *Gut Liver*. 2016;10(2):177-203.
72. Sebode M, Hartl J, Vergani D, Lohse AW, International Autoimmune Hepatitis G. Autoimmune hepatitis: From current knowledge and clinical practice to future research agenda. *Liver Int*. 2018;38(1):15-22.
73. Hennes EM, Zeniya M, Czaja AJ, Pares A, Dalekos GN, Krawitt EL, et al. Simplified criteria for the diagnosis of autoimmune hepatitis. *Hepatology*. 2008;48(1):169-76.
74. Kahaly GJ, Bartalena L, Hegedus L, Leenhardt L, Poppe K, Pearce SH. 2018 European Thyroid Association Guideline for the Management of Graves' Hyperthyroidism. *Eur Thyroid J*. 2018;7(4):167-86.
75. Persani L, Brabant G, Dattani M, Bonomi M, Feldt-Rasmussen U, Fliers E, et al. 2018 European Thyroid Association (ETA) Guidelines on the Diagnosis and Management of Central Hypothyroidism. *Eur Thyroid J*. 2018;7(5):225-37.
76. Bliddal S, Nielsen CH, Feldt-Rasmussen U. Recent advances in understanding autoimmune thyroid disease: the tallest tree in the forest of polyautoimmunity. *F1000Res*. 2017;6:1776.
77. Garmendia Madariaga A, Santos Palacios S, Guillen-Grima F, Galofre JC. The incidence and prevalence of thyroid dysfunction in Europe: a meta-analysis. *J Clin Endocrinol Metab*. 2014;99(3):923-31.
78. Tomer Y, Huber A. The etiology of autoimmune thyroid disease: a story of genes and environment. *J Autoimmun*. 2009;32(3-4):231-9.
79. Saravanan P, Dayan CM. Thyroid autoantibodies. *Endocrinol Metab Clin North Am*. 2001;30(2):315-37, viii.
80. Zophel K, Roggenbuck D, Schott M. Clinical review about TRAb assay's history. *Autoimmun Rev*. 2010;9(10):695-700.
81. Tozzoli R, Villalta D, Kodermaz G, Bagnasco M, Tonutti E, Bizzaro N. Autoantibody profiling of patients with autoimmune thyroid disease using a new multiplexed immunoassay method. *Clin Chem Lab Med*. 2006;44(7):837-42.
82. Ringel MD, Nabhan F. Approach to follow-up of the patient with differentiated thyroid cancer and positive anti-thyroglobulin antibodies. *J Clin Endocrinol Metab*. 2013;98(8):3104-10.

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Thermo Fisher Scientific – Phadia GmbH, Munzinger Str. 7, D-79111 Freiburg, Germany, Tel:
+49 761 47-805-0 © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are
the property of Thermo Fisher Scientific and its subsidiaries. Legal Manufacturer: Phadia AB,
Uppsala, Sweden 123808.AI.EU.EN.v1.21

ThermoFisher
S C I E N T I F I C